



Paper Chromatography Method - Sequencing Cards

Cut out the cards below and reorganise them into the correct sequence.

Lay the paper flat and leave it to dry.

Suspend the filter paper in a beaker so it isn't touching the sides or the bottom.

When the solvent has almost reached the top of the filter paper, remove it from the beaker.

Use a ruler to measure the distance the ink spot has travelled.

Place a spot of ink on the pencil line.

Cover the beaker and wait a few minutes.

Carefully add the solvent up to a level just below the ink spot, taking care not to splash the paper.

Draw a line with a pencil and ruler near to one end of the strip of filter paper.



Paper Chromatography Method Sequencing Cards - **Answer Sheet**

Draw a line with a pencil and ruler near to one end of the strip of filter paper.

Place a spot of ink on the pencil line.

Suspend the filter paper in a beaker so it isn't touching the sides or the bottom.

Carefully add the solvent up to a level just below the ink spot, taking care not to splash the paper.

Cover the beaker and wait a few minutes.

When the solvent has almost reached the top of the filter paper, remove it from the beaker.

Lay the paper flat and leave it to dry.

Use a ruler to measure the distance the ink spot has travelled.